



Anti-reverse wheel

A touch of extra safety for rehabilitation equipment

Anti-reverse wheel, as its name suggests, rotates only in one direction. It is an excellent choice for physical rehabilitation equipment and walking aids, among others. Anti-reverse is a safety feature to prevent the equipment from rolling backwards and stumbling the user which can lead to serious physical injuries.

The wheel is suitable for use with Tango forks and can be mounted on either a fixed or a directional lock fork.

The anti-reverse wheel is available in 125 mm size. Load capacity is similar to those of standard wheels. The wheel is available also with other tread hardnesses and materials as special products.



Technical information

Wheel diameter	125 mm
Fork material	Glass fibre reinforced polyamide
Rated load	100 kg
Fork options	Fixed and swivel with directional lock
Wheel material	Polypropylene centre and thermoelastic rubber tyre
Bearing of the wheel	Sealed ball bearing
Tread hardness	80 Shore A
Markings	Grey round sticker on the wheel centre Grey cover plate



Oy Mannerin Konepaja Ab, P.O.Box 3, 10901 Hanko, Finland
Tel. +358 19 222 001, Fax +358 19 248 2000
E-mail: manner@manner.fi, Internet: www.manner.fi

